

Raize Foot

The Raize Foot is lightweight and low profile microprocessor ankle/foot system. It's plantar/dorsi-flexion range decreases harmful forces on the residual limb and improves stability on slopes and slippery surfaces. By pressing a single button the wearer can control rollover resistance, lock or unlock the ankle, and manually adjust the heel height as needed.



Features and Benefits

Microprocessor controlled system

Easily adjustable heel height and lock (via control pad)

Reduces force on the residual limb

Accommodates various terrains

Automatically adjusts to allow maximum dorsiflexion on incline

Fillauer®



Raize Foot

The Raize foot is a versatile ankle/foot system with the ability to adapt to various terrains through its unique microprocessor system. The powerful system simulates natural ankle movement and reduces vertical shock, providing a comfortable walking gait. The easy to use control pad allows the user to modify plantar flexion, dorsiflexion range, and resistance whilst also providing complete control of the heel height.

The Raize foot is suitable for patients who partake in moderate everyday activity and exercise that requires prosthetic use beyond simple locomotion.

PROFILE INFORMATION

Activity Level	Mild to Moderate - (K2 and K3)		
Weight Limit	Up to 100 kg	Warranty	Foot 24 month limited Foot Shell 6 months
Size	24 - 30 cm	Foot Shell	Caucasian Brown Dark Brown
Heel Height Adjustment Dorsiflexion (up to 10°) Plantar flexion (up to 18°)	24 - 26 = 44 mm 27 - 28 = 51 mm 29 - 30 = 57 mm		

BUILD HEIGHT

Size	Build Height
24-28	89 mm
29-30	95 mm

Part Number Creation Raize Foot

Example 681-062-877

Product Code	Category Chart	Foot Size						
		24cm	25cm	26cm	27cm	28cm	29cm	30cm
682-062								
45 - 59 kg	869	872	876	880	884	917	N/A	
60 - 73 kg	870	873	877	881	885	918	921	
74 - 86 kg	871	874	878	882	886	919	922	
87 - 100 kg	**	875	879	883	916	920	923	

Part Number Creation Footshell

Example 45-26-13-03-3

Product Code	Foot Size	Color
Left = 45	24cm = 24	Caucasian = 03
Right = 46	25cm = 25	Brown = 09
	26cm = 26	Dark Brown = 13
	27cm = 27	
	28cm = 28	
	29cm = 29	
	30cm = 30	

