

Patient ID: _____

Order Number: _____

Clinician: _____

Facility: _____

Patient Sex: M F *

*Male patients will be provided with blue as a default and female patients will be provided with pink as default



Standard Order Specifications

Wide

JumpStart Cricket (4" - 8")

Outer Shell

- Polyethylene for sizes 4.00 – 8.00 in; blue or pink

JumpStart Leap Frog (4" - 8")

Outer Shell

- Polyethylene for sizes 4.00 – 8.00 in; blue or pink

Straps

- Instep & forefoot: riveted layover

JumpStart Bunny (4" - 9")

Outer Shell

- Polyethylene for sizes 4.00 – 5.75 in; blue or pink
- Co-poly for sizes 6.00 – 9.00 in; white (blue or pink straps)

Straps

- Instep & forefoot: riveted layover
- Proximal strap: elastic inside/outside strap

JumpStart Kangaroo (4" - 9")

Outer Shell

- Polyethylene for sizes 4.00 – 8.00 in; blue or pink
- Co-poly for sizes 8.25 – 9.00 in; white (blue or pink straps)

Inner Liner

- Height is same as foot length

Straps

- Instep: riveted d-ring
- Forefoot: riveted layover

Special Order Specifications

Narrow

8.25" - 9" JumpStart Cricket/Leapfrog

- Co poly for sizes 8.25 - 9.00 in; white

Adaptations

- Straps
- Padding
- Lining

Any requests not addressed on this form can be put in the Special Instructions.

Shoe Size (Euro)	Shoe Size (UK)	Fast Fit Sizing
15	0	4"
16	1	4.25"
17	1.5	4.50"
18	2.5	4.75"
19	3	5"
20	4	5.25"
21	5	5.50"
22	5.5	5.75"
23	6	6"
24	7	6.25"
25	8	6.5"
26	8.5	6.75"
27	9.5	7"
28	10	7.25"
29	11	7.50"
30	12	7.75"
31	12.5	8"

Size	Colour	Cricket®			Leap Frog®			Bunny®			Kangaroo®			
		Left	Right	Pair	Left	Right	Pair	Left	Right	Pair	Left	Right	Pair	
4.00 in	Blue													
10.2 cm	Pink													
4.25 in	Blue													
10.8 cm	Pink													
4.50 in	Blue													
11.4 cm	Pink													
4.75 in	Blue													
12.0 cm	Pink													
5.00 in	Blue													
12.7 cm	Pink													
5.25 in	Blue													
13.3 cm	Pink													
5.50 in	Blue													
14.0 cm	Pink													
5.75 in	Blue													
14.6 cm	Pink													
6.00 in	Blue													
15.3 cm	Pink													
6.25 in	Blue													
15.9 cm	Pink													
6.50 in	Blue													
16.5 cm	Pink													
6.75 in	Blue													
17.2 cm	Pink													
7.00 in	Blue													
17.8 cm	Pink													
7.25 in	Blue													
18.4 cm	Pink													
7.50 in	Blue													
19.0 cm	Pink													
7.75 in	Blue													
19.7 cm	Pink													
8.00 in	Blue													
20.3 cm	Pink													
8.25 in	Blue													
20.9cm	Pink													
8.50 in	Blue													
21.5 cm	Pink													
8.75 in	Blue													
22.2 cm	Pink													
9.00 in	Blue													
22.8cm	Pink													

Special Instructions