

Genudyn® CI NOVEL

Uniquely fastened for additional knee comfort

The Genudyn CI NOVEL brace is an aluminium knee orthosis with a unique, easily adjustable fastening system for control and dynamic stabilisation of the ACL and PCL.

The innovative twist-lock NOVEL system allows the strap position and tension to be permanently set, unlike most knee braces which rely on hook and loop fastenings that can shift with movement. This system is easily adjusted via tactile and audible check, even through clothing, allowing the user to release tension with a simple button press and remain consistently comfortable.



Features and benefits

- Extension adjustment: 0°, 10°, 20°, 30°
- Flexion adjustment: 0°, 10°, 20°, 30°, 45°, 60°, 90°
- Immobilisation possible at 0°, 10°, 20°, 30°
- Aluminium frame can be individually shaped, guaranteeing an accurate anatomical fit
- Ventral frame allows easy positioning without stepping in
- Incongruent joints support the knee's physiological movements
- Living Hinge provides compensation for frame motion
- Joint can be limited without the need for tools
- Functional knee guidance with stabilisation
- Additional strapping can be used to control and stabilise the posterior drawer
- Particularly suitable for active wearers, especially during sport

Indications for use

- Lesions, ruptures, and instability of the ACL, PCL, and lateral ligaments
- Injuries and reconstruction of the ACL, PCL, or both
- Instability of the anterior and posterior drawer
- Conservative/post-operative
- Combination injuries
- Meniscal injuries
- Patella fracture



Ordering information

Part number	Size	Circumference 15cm above the knee joint gap
Left knee		
OE-7785-XS-L	XS	32 - 38cm
OE-7785-S-L	S	38 - 44cm
OE-7785-M-L	M	44 - 50cm
OE-7785-L-L	L	50 - 56cm
OE-7785-XL-L	XL	56 - 63cm
OE-7785-XXL-L	XXL	63 - 70cm
Right knee		
OE-7785-XS-R	XS	32 - 38cm
OE-7785-S-R	S	38 - 44cm
OE-7785-M-R	M	44 - 50cm
OE-7785-L-R	L	50 - 56cm
OE-7785-XL-R	XL	56 - 63cm
OE-7785-XXL-R	XXL	63 - 70cm

